### Golden Horizon (Chengdu) Technology Co.,Ltd.

### **Technical Data Sheet**

Product Name: Chaga Extract

**Botanical source**: *Inonqqus obliquus* **Active ingredients**: Polysaccharides

Specification: 30%

Part used: Whole body

Appearance: Dark Brown Fine Powder

Test method: UV



#### **Health benefits:**

- 1) Chaga Mushroom Extract has been used to remedy in the treatment of stomach-intestinal tract diseases and as a palliative remedy for tumours of different location.
- 2) Chaga Mushroom Extract is used to cure skin diseases, especially in the case when they are combined with inflammatory diseases of the stomach-intestinal tract, liver and biliary colic.
- 3) Chaga Mushroom Extract improves health, preserves youth and increases longevity and boosts "Chi" or life energy.
- 4) It has been proved on numerous occasions, that using chaga products effectively helps to get rid of stones and sand in kidneys (nephrolith)

#### Application:

- 1) Chaga mushroom extract can be widely used in oral liquid,tablet,soft capsule,health care,pharmacy,beverage,cosmetics,food and other fields.
- 2) Chaga Extract Powder has been used to remedy in the treatment of stomach-intestinal tract diseases and as a palliative remedy for tumors of different location
- 3) Chaga Extract Powders used to cure skin diseases, especially in the case when they are combined with inflammatory diseases of the stomach-intestinal tract, liver and biliary colic.

# Golden Horizon (Chengdu) Technology Co.,Ltd.

## **Technical specification:**

ITEM	SPECIFICATION
Appearance	Dark-brown Fine Powder
Odor & taste	Characteristics
Particle size	100% through 80 meshes
Loss on drying	<b>≤7.0 %</b>
Ash	≤9.0%
Assay(Polysaccharides)	≥30%
As	≤1.0 mg/kg
Pb	≤2.0 mg/kg
Hg	≤0.2 mg/kg
Cd	≤1.0 mg/kg
Microbiological Analysis	
Total of Plate Count	<10000CFU/g
Yeast and Mold	<100CFU/g
E.Coli	Negative
Staphylococcus	Negative

## **Packaging Details:**

Package: 25kgs/fiber drum, inner with double-layer plastic bags.

Storage: Store in cool, dry place, and keep away from direct strong light.

**END OF DOCUMENT**